



ANTI HOLOGRAM



EAN code: 5907588408714 item no.: 0035121 cap. 250ml

EAN code: 5907588408738 item no.: 0035122 cap. 500ml

+ PROPERTIES

ANTI HOLOGRAM is a product for the most demanding users. It complements the BOLL 1000 toothpaste and ensures that the result is of the highest possible quality. It is a delicately abrasive paste that eliminates micro scratches and possible holograms created in previous polishing works. After applying the BOLL B1000 paste and the combination with a properly selected sponge, you can achieve the so-called "mirrors effect".

••• COLOUR white**↪** METHOD OF USE

Apply a small amount of paste on a properly selected polishing sponge attached to the polishing machine. Then, using a sponge, spread the paste on the part and start polishing, starting from the lowest speed. Gradually increase the speed of the polishing machine, making sure not to exceed 2 000 rpm. The best effect can be obtained after prior application of B1000 paste.

☰ NOTES

Do not use on hot elements!

↓ APPLICATION

Mix thoroughly



Can be used by hand and machine



Maximum 2000 rpm



Wipe with microfiber



Without silicone



PHYSICAL PROPERTIES

density at 20°C: 1,01 g/cm³
flash point: > 93°C



STORAGE CONDITIONS



Protect the product from excessive heat and cold. Avoid direct exposure to sunlight. Store in a dry place.



SHELF LIFE

24 months from the date of production placed on the package.



SAFETY

See Safety Data Sheet.



OTHER INFORMATION

All technical data are approximate values. We advise you to test the material to ensure suitability for your specific application. The producer reserves the right to improve the product and change the technical conditions with the possibility of making changes inside the specification.